

# New additions to the Labware family - completing our labware portfolio

Introducing the new Sperm collection container and 5 mL Sample tube – designed, developed and manufactured specifically for use in IVF.

## New Sperm collection container

The new Sperm collection container has been designed with the user in mind. With its oval shape and snap-on lid it enables easy and safe sampling. The first sampling cup specifically CE-marked for IVF.

## New 5 mL Sample tube

A small tube for your IVF lab with a two-position cap enabling conditioning of media. Sample tube

5 mL is a high quality polystyrene tube, CE-marked for IVF.

## Developed, designed, manufactured and certified for IVF

Vitrolife's high quality labware products are specifically designed for IVF and together they will enable you to create a secure, protective and supportive environment for gametes and embryos.

## Safe environment from aspiration to transfer

Vitrolife Labware gives you the possibility to use a high quality, extensively tested environment for your gametes and embryos through all the steps of an IVF procedure.



- Guaranteed safe environment for embryos and gametes
- CE-marked for IVF

# Safe, efficient and certified



## MEA-tested

Multiple endpoints, including 1-cell, expanded blastocysts within 96h  $\geq$  80% and cell count.



## HSSA-tested

Motility  $\geq$  80%.  
Products for sperm handling.



## 100% control

All stages of production are controlled, from the raw materials to the final product. Extensive sampling of each LOT.



## Sterile

Sterility assurance level (irradiation)  $10^{-6}$ .  
Less than 0.25 endotoxin units/device.



## Sperm collection container REF 16102

### Designed for you

- Supporting procedure and efficient workflow

### Meeting regulatory requirements

- First sperm sampling cup CE-marked for IVF

### Features

- 100 individually packed per case
- Unique design
- Conical shaped bottom for easy access with Pasteur pipette

- Tamper evident seal
- Snap on lid
- Can be used for IUI, sperm donation, IVF and ICSI
- Material: polypropylene



## 5 mL Sample tube REF 16103

### Designed for you

- Supporting procedure and efficient workflow

### Meeting regulatory requirements

- First 5 ml Sample tube CE-marked for IVF

### Features

- Manufactured, certified and tested for IVF
- 400 per case and 25 pcs/tray
- Two-position cap (PE) to enable venting
- Material: polystyrene

Contact your sales representative  
or order online at [www.vitrolife.com](http://www.vitrolife.com).

Note: All products shown on this flyer might not be available in all markets.

Vitrolife • Box 9080 • SE-400 92 Gothenburg • Sweden • Tel +46 31 721 80 00 • [order@vitrolife.com](mailto:order@vitrolife.com) • [www.vitrolife.com](http://www.vitrolife.com)